

## Product Datasheet

### Human IKBalpha protein (orb633416)

<b>Description</b>	Recombinant human IKBalpha protein
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB
<b>Preservatives</b>	Lyophilized from a 0.2 µm filtered solution in PBS with 5% trehalose, pH7.4
<b>Storage</b>	-20°C for 12 months as lyophilized;2-8°C for 1 month under sterile conditions after reconstitution
<b>Tag</b>	N-terminal GST tag
<b>Note</b>	For research use only
<b>Application notes</b>	Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 µg/ml in sterile distilled water by gentle pipetting 2-3 times, don't vortex
<b>Formula</b>	Lyophilized from a 0.2 µm filtered solution in PBS with 5% trehalose, pH <sub>7</sub> .4
<b>Purity</b>	Greater than 95% by SDS-PAGE gel analyses
<b>Source</b>	E.Coli
<b>MW</b>	57.4 KDa
<b>Uniprot ID</b>	<b>P25963</b>
<b>Expiration Date</b>	6 months from date of receipt.