

Product Datasheet

4E-BP1 (phospho-Ser65) (orb636064)

Description	Mouse monoclonal antibody to 4E-BP1 (phospho-Ser65)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IHC
Immunogen	Synthesized peptide of human Phospho-4E-BP1 (Ser65).
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	12.6kDa
Entrez	1978
Dilution Range	IHC: 1/200 - 1/1000; ICC: 1/100 - 1/500; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.