

Product Datasheet

Human Iba1 / AIF-1 Protein (orb636715)

Description	Human Iba1, His Tag is expressed from E. coli cells. It contains AA Met 1 - Pro 147 (Accession # P55...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Iba1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 18.6 kDa. The protein migrates as 18-20 kDa under reducing (R) condition (SDS-PAGE).
Purity	95%
Source	Human Iba1, His Tag (orb636715) is expressed from E. coli cells. It contains AA Met 1 - Pro 147 (Accession # P55008-1).
MW	18.6 kDa
Uniprot ID	P55008
NCBI	NP_001614.3
Expiration Date	6 months from date of receipt.