

## Product Datasheet

### Human LILRB4 / CD85k / ILT3 Protein (orb636719)

<b>Description</b>	Human LILRB4, His Tag is expressed from human 293 cells (HEK293). It contains AA Gln 22 - Glu 259 (A...
<b>Reactivity</b>	Human
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	LILRB4
<b>Preservatives</b>	PBS, pH7.4
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	C-10×His
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 28.2 kDa. The protein migrates as 32-43 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	AAH26309.1
<b>Purity</b>	90%
<b>Source</b>	Human LILRB4, His Tag (orb636719) is expressed from human 293 cells (HEK293). It contains AA Gln 22 - Glu 259 (Accession # AAH26309.1).
<b>MW</b>	28.2 kDa
<b>NCBI</b>	<a href="#">AAH26309.1</a>
<b>Expiration Date</b>	6 months from date of receipt.