

Product Datasheet

Human PSCA Protein (orb636720)

Description	Human PSCA, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 12 - Ser 86 (Acce...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PSCA
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 10.4 kDa. The protein migrates as 18-25 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_005663.2
Purity	90%
Source	Human PSCA, His Tag (orb636720) is expressed from human 293 cells (HEK293). It contains AA Leu 12 - Ser 86 (Accession # O43653-1).
MW	10.4 kDa
Uniprot ID	O43653
NCBI	NP_005663.2
Expiration Date	6 months from date of receipt.