

## Product Datasheet

### Human ErbB3 / Her3 Protein (orb636723)

<b>Description</b>	Human ErbB3 Protein, Twin-Strep Tag is expressed from human 293 cells (HEK293). It contains AA Ser 2...
<b>Reactivity</b>	Human
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	ErbB3
<b>Preservatives</b>	PBS, pH7.4
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	C-Twin-Strep
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a twin strep tag at the C-terminus. The protein has a calculated MW of 71.7 kDa. The protein migrates as 85-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	NP_001973.2
<b>Purity</b>	95%
<b>Source</b>	Human ErbB3, Strep Tag (orb636723) is expressed from human 293 cells (HEK293). It contains AA Ser 20 - Thr 643 (Accession # NP_001973.2).
<b>MW</b>	71.7 kDa
<b>Uniprot ID</b>	<a href="#">P21860</a>
<b>NCBI</b>	<a href="#">NP_001973.2</a>
<b>Expiration Date</b>	6 months from date of receipt.