
Product Datasheet

Human ErbB3 / Her3 Protein (orb636724)

Description

Human ErbB3, Fc Tag is expressed from

Ser 20 - Thr 643 (Acc...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

ErbB3

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 94.8 kDa. The protein migrates as 130-170 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_001973.2

Purity

95%

Source

Human ErbB3, Fc Tag (orb636724) is expressed from human 293 cells (HEK293). It contains AA Ser 20 - Thr 643 (Accession # NP_001973.2).

MW

94.8 kDa

Uniprot ID

[P21860](#)

NCBI

[NP_001973.2](#)

Expiration Date

6 months from date of receipt.