

## Product Datasheet

### IPA-3 (orb636833)

|                          |                                  |
|--------------------------|----------------------------------|
| <b>Description</b>       | IPA-3                            |
| <b>Target</b>            | Kinase                           |
| <b>Form/Appearance</b>   | powder                           |
| <b>Storage</b>           | -20°C                            |
| <b>Note</b>              | For research use only            |
| <b>Formula</b>           | $C_{20}H_{14}O_2S_2$             |
| <b>Purity</b>            | > 98%                            |
| <b>MW</b>                | 350.45                           |
| <b>CAS Number</b>        | 42521-82-4                       |
| <b>Solubility (25°C)</b> | Soluble in DMSO (up to 25 mg/ml) |
| <b>Expiration Date</b>   | 12 months from date of receipt.  |