

Product Datasheet

AKT1/2/3 (Phospho-Tyr315/Tyr316/Tyr312) antibody (orb636869)

Description	Rabbit polyclonal antibody to AKT1/2/3 (Phospho-Tyr315/Tyr316/Tyr312)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1: 100-1:500, ICC: 1: 100-1:500
Expiration Date	12 months from date of receipt.