


Product Datasheet

Torc2 (phospho-Ser171) antibody (orb6371)

Description	Torc2 (phospho-Ser171) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human CREB regulated transcription coactivator 2 around the phosphorylation site of Ser171 (TS(p-S)DS)
Target	Phospho-TORC2 (Ser171)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	73 kDa
Uniprot ID	Q53ET0
Dilution Range	WB=1:500-2000, Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.



Western blot analysis of human rectal ca...