

## Product Datasheet

### TM158 antibody (orb637511)

<b>Description</b>	Rabbit polyclonal antibody to TM158
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IF, WB
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	WB: 1:500-1:2000, IF: 1: 100-1:500, ICC: 1: 100-1:500
<b>Expiration Date</b>	12 months from date of receipt.