

Product Datasheet

ARHGAP10 antibody (orb637809)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to ARHGAP10 |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | ICC, IF, WB |
| Note | For research use only |
| Clonality | Polyclonal |
| Dilution Range | WB: 1:500-1:2000, IF: 1: 100-1:500, ICC: 1: 100-1:500 |
| Expiration Date | 12 months from date of receipt. |