

## Product Datasheet

# Golgin subfamily A member 4 antibody (orb638252)

<b>Description</b>	Rabbit polyclonal antibody to Golgin subfamily A member 4
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, ICC, IF, IHC-P, WB
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	IF: 1: 100-1:500, ICC: 1: 100-1:500, WB: 1:500-1:2000, IHC-P: 1:50-1:200, ELISA: 1:20000-1:40000
<b>Expiration Date</b>	12 months from date of receipt.