

Product Datasheet

B7-H6 antibody (orb638494)

Description	Rabbit polyclonal antibody to B7-H6
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IF: 1: 100-1:500, ICC: 1: 100-1:500, WB: 1:500-1:2000, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.