

Product Datasheet

Flag-Tag antibody (orb638733)

Description	Rabbit polyclonal antibody to Flag-Tag
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IP, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1: 1000-1:30000, IF: 1:200, ICC: 1:200, IP 1:200, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.