
Product Datasheet

Human BTLA (31-150) Protein (orb638745)

Description

Human BTLA (31-150), Mouse IgG2a Fc

(HEK293). It contains AA L...

Reactivity

Human

Endotoxins

0.1 EU per µg

Conjugation

Unconjugated

Target

BTLA

Preservatives

50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

C-mFc(mlgG2a)

Note

For research use only

Application notes

This protein carries a mouse IgG2a Fc tag at the C-terminus. The protein has a calculated MW of 40.6 kDa. The protein migrates as 50-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

AAP44003.1

Purity

95%

Source

Human BTLA (31-150), Mouse IgG2a Fc Tag (orb638745) is expressed from human 293 cells (HEK293). It contains AA Lys 31 - Ser 150 (Accession # AAP44003.1).

MW

40.6 kDa

Uniprot ID
[Q7Z6A9](#)
NCBI
[AAP44003.1](#)
Expiration Date

6 months from date of receipt.