

---

## Product Datasheet

### Mouse Properdin Protein (orb638749)

**Description**

Mouse Properdin, His Tag is expressed from

Ser 23 - Pro 464...

**Reactivity**

Mouse

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

Properdin

**Preservatives**

PBS, pH7.4

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

**Tag**

C-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 49.8 kDa. The protein migrates as 55-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

95%

**Source**

Mouse Properdin, His Tag (orb638749) is expressed from human 293 cells (HEK293). It contains AA Ser 23 - Pro 464 (Accession # P11680-1).

**MW**

49.8 kDa

**Uniprot ID**

[P11680](#)

**NCBI**

[NP\\_032849.2](#)

**Expiration Date**

6 months from date of receipt.