
Product Datasheet

Human PSMA / FOLH1 Protein (orb638750)

Description

Human PSMA, Fc Tag is expressed from

Lys 44 - Ala 750 (Acce...

Reactivity

Human

Endotoxins

0.01 EU per µg

Conjugation

Unconjugated

Target

PSMA

Preservatives

50 mM Tris, 100 mM Glycine, 150 mM NaCl, pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

N-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the N-terminus. The protein has a calculated MW of 106.0 kDa. The protein migrates as 110-140 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_058050.3

Purity

95%

Source

Human PSMA, Fc Tag (orb638750) is expressed from human 293 cells (HEK293). It contains AA Lys 44 - Ala 750 (Accession # Q04609-1).

MW

106.0 kDa

Uniprot ID

[Q04609](#)

NCBI

[NP_058050.3](#)

Expiration Date

6 months from date of receipt.