

---

## Product Datasheet

### Cynomolgus TSLP Protein (orb638752)

**Description**

Cynomolgus TSLP, His Tag is expressed

AA Tyr 29 - Gln 159...

**Reactivity**

Monkey

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

TSLP

**Preservatives**

PBS, pH7.4

**Form/Appearance**

Powder

**Storage**

-20°C

**Tag**

C-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 17.1 kDa. The protein migrates as 8-10 kDa, 15-20 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Protein Sequence**

XP\_005557555.1

**Purity**

90%

**Source**

Cynomolgus TSLP, His Tag (orb638752) is expressed from human 293 cells (HEK293). It contains AA Tyr 29 - Gln 159 (Accession # A0A7N9CAT7-1).

**MW**

17.1 kDa

**Uniprot ID**

[A0A7N9CAT7](#)

**NCBI**

[XP\\_005557555.1](#)

**Expiration Date**

6 months from date of receipt.