



Product Datasheet

SARS N (nucleocapsid) antibody (orb638865)

Description	SARS N (nucleocapsid) antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant SARS N (nucleocapsid) (1-422/422aa)
Target	SARS-CoV Nucleocapsid
Preservatives	0.01M PBS(pH7.4) with 0.1% Proclin300
Form/Appearance	Liquid
Concentration	≥1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	46 kDa
Uniprot ID	P59594
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1:2019-nCoV Nucleocaps...



Western blot: Lane 1:SARS-CoV-2 N Protein...