


Product Datasheet

SARS-CoV-2 N protein antibody (orb638873)

Description	SARS-CoV-2 N protein antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid Protei (full length)
Target	SARS-CoV-2 (2019-nCoV) Nucleocapsid
Preservatives	0.01M PBS(pH7.4) with 0.1% Proclin300
Form/Appearance	Liquid
Concentration	≥1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	46 kDa
Uniprot ID	PODTC2
Dilution Range	WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: SARS-CoV-2 N protein...