

Product Datasheet

CCR1 antibody (orb639028)

Description	Rabbit polyclonal antibody to CCR1
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human CCR1(NP_001286.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41kDa
Uniprot ID	P32246
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.