

Product Datasheet

TAF2 antibody (orb639068)

Description	Rabbit polyclonal antibody to TAF2
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TAF2(NP_003175.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6P1X5
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.