

Product Datasheet

RDH11 antibody (orb639135)

Description	Rabbit polyclonal antibody to RDH11
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human RDH11(NP_057110.3).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35kDa
Uniprot ID	Q8TC12
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.