

Product Datasheet

ASK1 (phospho-Ser1033) antibody (orb6392)

Description	ASK1 (phospho-Ser1033) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human ASK1 around the phosphorylation site of Ser1033 (PP(p-S)PE)
Target	phospho-ASK1 (Ser1033)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	155 kDa
Uniprot ID	Q99683
Expiration Date	12 months from date of receipt.