

Product Datasheet

MMP17 antibody (orb639335)

Description	Rabbit polyclonal antibody to MMP17
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human MMP17(NP_057239.4).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	75kDa
Uniprot ID	Q9ULZ9
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.