

Product Datasheet

MAP2K7 (Phospho-S271/T275) antibody (orb639379)

Description	Rabbit polyclonal antibody to MAP2K7 (Phospho-S271/T275)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S271 & T275 of human MEK7(NP_660186.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47kDa/48kDa/49kDa/51kDa
Uniprot ID	O14733
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.