

Product Datasheet

ATP1A1 (Phospho-S16) antibody (orb639390)

Description	Rabbit polyclonal antibody to ATP1A1 (Phospho-S16)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S16 of human Na ⁺ /K ⁺ -ATPase(NP_000692.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	112kDa
Uniprot ID	P05023
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.