

Product Datasheet

AKT1 (Phospho-S473) & AKT2-(Phospho-S474/AKT3-S472) antibody (orb639392)

Description	Rabbit polyclonal antibody to AKT1 (Phospho-S473) & AKT2-(Phospho-S474/AKT3-S472)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A phospho synthetic peptide of human AKT1/AKT2/AKT3.
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	60kDa
Uniprot ID	P31749
Entrez	20720810000
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.