

Product Datasheet

Bcl (Phospho-2 (Phospho-S87) antibody (orb639395)

Description	Rabbit polyclonal antibody to Bcl (Phospho-2 (Phospho-S87)
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S87 of human Bcl-2(NP_000624.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28kDa
Uniprot ID	P10415
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.