

Product Datasheet

FOXO1 (Phospho-S319) antibody (orb639411)

Description	Rabbit polyclonal antibody to FOXO1 (Phospho-S319)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic phosphorylated peptide around S319 of human FOXO1(NP_002006.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90kDa
Uniprot ID	Q12778
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.