

Product Datasheet

PTK2 (Phospho-Y576/577) antibody (orb639416)

Description	Rabbit polyclonal antibody to PTK2 (Phospho-Y576/577)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y925 of human FAK(NP_722560.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39kDa/48kDa/63kDa/99kDa/114kDa/119kDa/120kDa
Uniprot ID	Q05397
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.