

Product Datasheet

RPS6KB1 (Phospho-Thr421) antibody (orb639423)

Description	Rabbit polyclonal antibody to RPS6KB1 (Phospho-Thr421)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T421 of human RPS6KB1(NP_001258971.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	51kDa/52kDa/56kDa/59kDa
Uniprot ID	P23443
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.