

Product Datasheet

STK11 (Phospho-S428) antibody (orb639426)

Description	Rabbit polyclonal antibody to STK11 (Phospho-S428)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding S428 of human STK11.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	65kDa
Uniprot ID	Q15831
Entrez	6794
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.