

Product Datasheet

p53 (Phospho-S37) antibody (orb639428)

Description	Rabbit polyclonal antibody to p53 (Phospho-S37)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic phosphorylated peptide around S37 of human p53(NP_000537.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	53KDa
Uniprot ID	P04637
Dilution Range	IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.