

Product Datasheet

EGFR (Phospho-Y869) antibody (orb639431)

Description	Rabbit polyclonal antibody to EGFR (Phospho-Y869)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y869 of human EGFR(NP_005219.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	134KDa
Uniprot ID	P00533
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.