

Product Datasheet

PAK1 (Phospho-S199/S204) antibody (orb639442)

Description	Rabbit polyclonal antibody to PAK1 (Phospho-S199/S204)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding S199/S204 of human PAK1.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q13153
Entrez	5058
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.