

Product Datasheet

Tau (Phospho-Thr212/Ser214) antibody (orb639443)

Description	Rabbit polyclonal antibody to Tau (Phospho-Thr212/Ser214)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding Thr212/Ser214 of human Tau.
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	50kDa
Uniprot ID	P10636
Entrez	4137
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.