

Product Datasheet

Pan Cytokeratin Cocktail (Acidic + Basic) antibody (FITC) (orb639460)

Description

Mouse monoclonal antibody against

Basic) conjugated to FITC

Species/Host

Mouse

Reactivity

Human

Conjugation

FITC

Tested

FACS, IF

Applications
Immunogen

Human epidermal keratin was used as the immunogen for this Pan Cytokeratin antibody.

Preservatives

500 ul at 0.1 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the Pan Cytokeratin antibody at 2-8°C, protected from light.

Note

For research use only

Application notes

Optimal dilution of the Pan Cytokeratin antibody should be determined by the researcher.

Formula

500 ul at 0.1 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Hazard
Information

This Pan Cytokeratin antibody is available for research use only.

Dilution Range

Flow cytometry: 5ul per test per one million cells in 0.1ml or 5ul per 100ul of whole blood, Immunofluorescence: 1:50-1:100

Expiration Date

12 months from date of receipt.