

Product Datasheet

SPTBN2 Antibody / Spectrin beta III (orb639749)


Description

Spectrin is an actin binding protein that is a major


	Spec...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	A recombinant human partial protein (amino acids 356-475) was used as the immunogen for the recombinant Spectrin beta III antibody.
Preservatives	0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the recombinant Spectrin beta III antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the antibody should be determined by the researcher.
Formula	0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Recombinant
Purity	Protein G affinity chromatography
Uniprot ID	O15020
Hazard Information	This recombinant Spectrin beta III antibody is available for research use only.
Dilution Range	ELISA (order BSA-free format for coating), Western blot: 1-2ug/ml
Expiration Date	12 months from date of receipt.



Western blot testing of human HeLa cell ...



Analysis of HuProt (TM) microarray conta...



SDS-PAGE analysis of purified, BSA-free ...