
Product Datasheet

Calcineurin Antibody / PPP3CA (orb639860)

Description

Calcineurin is an enzyme that

residues in proteins. It is a he...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

WB

Immunogen

A recombinant full-length Calcineurin A protein was used as the immunogen for the PPP3CA antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the PPP3CA antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

[Q08209](#)

Hazard Information

This PPP3CA antibody is available for research use only.

Dilution Range

Western blot: 1-2ug/ml

Expiration Date

12 months from date of receipt.