


Product Datasheet

Calcineurin B Antibody / PPP3R1 (orb639878)

| | |
|----------------------------|--|
| Description | Calcineurin is an enzyme that dephosphorylates serine and threonine residues in proteins. It is a he... |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IHC-P |
| Immunogen | A recombinant full-length human Calcineurin B protein was used as the immunogen for the PPP3R1 antibody. |
| Preservatives | 0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide |
| Storage | Store the PPP3R1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide). |
| Note | For research use only |
| Application notes | Optimal dilution of the antibody should be determined by the researcher. |
| Formula | 0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide |
| Isotype | Mouse IgG1, kappa |
| Clonality | Monoclonal |
| Purity | Protein G affinity chromatography |
| Uniprot ID | P63098 |
| Hazard Information | This PPP3R1 antibody is available for research use only. |
| Dilution Range | Immunohistochemistry (FFPE): 1-2ug/ml |
| Expiration Date | 12 months from date of receipt. |



IHC staining of FFPE human prostate carc...