

---

## Product Datasheet

### Granulocyte-Macrophage CSF Antibody / CSF2 (orb639958)

**Description**

Granulocyte/macrophage - Colony-stimulating



produc...

**Species/Host**

Rat

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested Applications**

IHC-P

IHC staining of FFPE human spleen with G...

**Immunogen**

Recombinant human protein was used as the immunogen for this GMCSF antibody.

**Preservatives**

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

**Storage**

Store the GMCSF antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

**Note**

For research use only

**Application notes**

Optimal dilution of the antibody should be determined by the researcher.

IHC staining of FFPE human spleen with G...

**Formula**

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

**Isotype**

Mouse IgG2b, kappa

**Clonality**

Monoclonal

**Purity**

Protein G affinity chromatography

**Uniprot ID**
**P04141**
**Hazard Information**

This GMCSF antibody is available for research use only.

IHC staining of FFPE human spleen with G...

**Dilution Range**

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

**Expiration Date**

12 months from date of receipt.