

Product Datasheet

EBAG9 Antibody / Estrogen Receptor Binding Site Associated Antigen 9 (orb640025)

Description

EBAG9, also known as RCAS1, is an estrogen-


Species/Host

RCAS1 protein...

Mouse

Reactivity


Human

Conjugation

Unconjugated

Tested Applications

IF, IHC-P



Immunofluorescent staining of permeabilized cells...

Immunogen


A recombinant human full-length EBAG9 protein was used as the immunogen for this EBAG9 antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the EBAG9 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).



IHC staining of FFPE human tonsil tissue...

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype


Mouse IgG2c, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography



SDS-PAGE analysis of purified, BSA-free ...

Uniprot ID

O00559

Hazard Information

This EBAG9 antibody is available for research use only.

Dilution Range

Immunofluorescence: 1-2ug/ml, Immunohistochemistry (FFPE): 0.5-1ug/ml

Expiration Date

12 months from date of receipt.