



Product Datasheet

CD44v4/5 Antibody (orb640056)


Description	This antibody recognizes an epitope encoded by exon v9 on the variant portion of human CD44. The CD4...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IHC-P, WB
Immunogen	A recombinant human full-length CD44v3-10 protein was used as the immunogen for this CD44v4/5 antibody.
Preservatives	0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the CD44v4/5 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the antibody should be determined by the researcher.
Formula	0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P16070
Hazard Information	This CD44v4/5 antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells in 0.1ml,Western blot: 1-2ug/ml,Immunohistochemistry (FFPE): 1-2ug/ml
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human cervix with C...



IHC staining of FFPE human cervix with C...



SDS-PAGE analysis of purified, BSA-free ...