



Product Datasheet

NKX3.1 Antibody (orb640060)


Description	NKX3.1 is a prostate specific gene encoding a transcription factor that plays an important role in n...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	A recombinant human partial protein (amino acids 92-224) was used as the immunogen for the NKX3.1 antibody.
Preservatives	0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the NKX3.1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the antibody should be determined by the researcher.
Formula	0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	Q99801
Hazard Information	This NKX3.1 antibody is available for research use only.
Dilution Range	ELISA (order BSA-free format for coating), Immunohistochemistry (FFPE): 1-2ug/ml
Expiration Date	12 months from date of receipt.



SDS-PAGE analysis of purified, BSA-free ...



IHC staining of FFPE human prostate carc...



IHC staining of FFPE human prostate carc...