

Product Datasheet

Enterobacter aerogenes Antibody (orb640086)

Description	Enterobacter aerogenes is a Gram-negative rod-shaped microorganism from the Enterobacteriaceae family...
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	A crude sonicate of Enterobacter aerogenes was used as the immunogen for the Enterobacter aerogenes antibody.
Preservatives	0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the Enterobacter aerogenes antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the antibody should be determined by the researcher.
Formula	0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Hazard Information	This Enterobacter aerogenes antibody is available for research use only.
Dilution Range	ELISA (order BSA-free format for coating), Immunofluorescence: 1-2ug/ml
Expiration Date	12 months from date of receipt.