


Product Datasheet

Pan Cytokeratin Antibody Cocktail (Acidic + Basic) (orb640118)

Description

Anti-cytokeratin clone PCK/3150 demonstrates a



normal t...

Species/Host

Mouse

Reactivity

Human


Conjugation

Unconjugated

Tested

Applications

FACS, IF, IHC-P




Immunofluorescent
staining of
permeabili...

Immunogen

A crude cytokeratin extract prepared from RT-4 and MCF-7 cells was used as the immunogen for this Pan Cytokeratin antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide



Flow cytometry
testing of
permeabilized ...

Storage

Store the Pan Cytokeratin antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2a, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Hazard

Information


This Pan Cytokeratin antibody is available for research use only.

Dilution Range

Flow cytometry: 1-2ug/million cells in 0.1ml, Immunofluorescence: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml

Expiration Date

12 months from date of receipt.



IHC staining of
FFPE human
prostate with...