
Product Datasheet

Human TROP-2 / TACSTD2 Protein (orb640175)

Description

Human TROP-2, Fc Tag is expressed

contains AA Gln 31 - Thr 274 (Ac...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

TROP-2

PreservativesTris with Glycine, Arginine and NaCl,
pH7.5**Form/Appearance**

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 54.3 kDa. The protein migrates as 60-70 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_002344.2

Purity

95%

Source

Human TROP-2, Fc Tag (orb640175) is expressed from human 293 cells (HEK293). It contains AA Gln 31 - Thr 274 (Accession # P09758-1).

MW

54.3 kDa

Uniprot ID[P09758](#)**NCBI**[NP_002344.2](#)**Expiration Date**

6 months from date of receipt.