

Product Datasheet

DDDDK-tag (1A8) antibody (orb640818)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to DDDDK-tag (1A8) |
| Species/Host | Mouse |
| Reactivity | All |
| Conjugation | Unconjugated |
| Tested Applications | IF, IP, WB |
| Immunogen | Synthetic peptide DDDDK-tag. |
| Preservatives | 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2 |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Monoclonal |
| Dilution Range | WB: 1:5000-1:50000IP: 1:100-1:1000, IF: 1:100-1:1000 |
| Expiration Date | 12 months from date of receipt. |